

# A Novel Mixing Device For the Reproducible Manufacture of Non-Viral Gene Therapy Formulations

Lee A. Davies<sup>1,4</sup>, Graciela A. Nunez-Alonso<sup>1,4</sup>, Henry L. Hebel<sup>2</sup>, Ron K. Scheule<sup>3</sup>, Seng H. Cheng<sup>3</sup>, Deborah R. Gill<sup>1,4</sup> and Stephen C. Hyde<sup>1,4</sup>.

1. Gene Medicine Group, Nuffield Department of Clinical Laboratory Sciences, University of Oxford, John Radcliffe Hospital, Oxford, UK.

2. VGXI, Houston, TX, USA.

3. Genzyme Corporation, Framingham, MA, USA.

4. The UK Cystic Fibrosis Gene Therapy Consortium, Edinburgh/London/Oxford.

lee.davies@ndcls.ox.ac.uk

http://users.ox.ac.uk/~genemed

http://www.cfgenetherapy.org.uk

Poster download available

## INTRODUCTION

- ▶ Production of non-viral gene therapy formulations requires the controlled mixing of GTAs with plasmid DNA (pDNA)
- ▶ Variability in mixing procedures between laboratories and individuals can greatly impact formulation efficacy
- ▶ Large-scale production of non-viral gene therapy formulations suitable for clinical applications will require a reliable and reproducible mixing methodology
- ▶ We have developed and tested a novel pneumatic mixing device for the standardised production of non-viral gene therapy formulations

## RESULTS




Pneumatic piston  
Plunger  
Dual lumen syringe

Dimensions: w : d : h = 20 x 12 x 45cm


### The LMD2 pneumatic mixer device

- ▶ Fully adjustable pneumatic syringe-driver powered via 50 psi compressed air source
- ▶ Incorporates dual lumen syringe containing pDNA and GTA in separate compartments
- ▶ Activation drives both reagents through 8-element HDPE static mixer to ensure adequate mixing
- ▶ Fully scaleable for large volume clinical mixing



### Electron micrographs of experimental formulations

- ▶ Mixing studies were performed using two important non-viral GTAs that show promise for aerosol gene delivery
- ▶ GL67A (GL67:DOPE:DMPE-PEG<sub>5000</sub>) is a cationic lipid that forms complex multi-lamellar lipoplexes with pDNA
- ▶ 25 kDa branched polyethylenimine (PEI) is a highly cationic polymer that forms small discrete polyplexes with pDNA



Data represents mean ±SEM for n=3 separate measurements

### Reproducible mixing using the LMD2


- ▶ LMD2 mixer was loaded with GL67A and the luciferase expression plasmid pCIKLux (5.6 kb)
- ▶ Complexes were prepared using a range of applied mixing rates from 0.5 ml/s to 10 ml/s
- ▶ Vertical displacement of LMD2 plunger during activation measured using high-speed video analysis
- ▶ Linear displacement indicates constant and reproducible mixing conditions at all rates tested

		Mixing rate (ml/s)					
		Hand	0.5	1.0	2.5	5.0	10.0
GL67A*	Size (nm)	265 ± 18	287 ± 13	278 ± 13	285 ± 19	306 ± 16	315 ± 22
	Charge (mV)	+3.8 ± 0.2	+3.5 ± 0.3	+3.6 ± 0.4	+3.4 ± 0.3	+4.1 ± 0.5	+4.1 ± 0.4
PEI**	Size (nm)	75 ± 4	73 ± 2	70 ± 1	73 ± 1	80 ± 1	95 ± 2
	Charge (mV)	+47.3 ± 0.7	+47.1 ± 0.9	+48.3 ± 0.4	+48.8 ± 0.3	+49.9 ± 0.4	+48.4 ± 0.8

Data represents mean ±SEM for n=3 separate measurements.  
\*GL67A : pDNA ratio 0.6mM : 0.8mM. \*\*PEI : pDNA at N : P ratio of 10 : 1.


### Physical characteristics of generated complexes

- ▶ pCIKLux/GL67 and pCIKLux/PEI complexes prepared over a range of mixing rates from 0.5 ml/s to 10 ml/s or by hand
- ▶ Physical characteristics measured using dynamic light scattering (Zetasizer Nano, Malvern Instruments, UK)
- ▶ pCIKLux/GL67 complexes demonstrate equivalent size/charge characteristics to complexes prepared by hand
- ▶ pCIKLux/PEI complexes increase in size at higher mixing rates due to formation of aggregates



### Absence of shear degradation using the LMD2

- ▶ Plasmid DNA is highly sensitive to shear degradation during manufacturing processes
- ▶ pCIKLux integrity analysed by gel electrophoresis after complexation with GL67A or PEI using the LMD2
- ▶ No evidence of increased shear related damage following complex preparation using the LMD2 compared to standard hand mixing



Data represents mean ±SEM for n=3 separate experiments.  
GL67A : pDNA ratio 6mM : 8mM. PEI : pDNA at N : P ratio of 10 : 1.

### In vivo performance of LMD2 formulations

- ▶ 10 ml pCIKLux/PEI or pCIKLux/GL67 complexes generated using the LMD2 at a mixing rate of 2.5 ml/s
- ▶ Aerosolised to BALB/c mice (n=6) using whole body aerosol exposure chamber and pneumatic nebuliser
- ▶ Luciferase expression in lung tissues analysed 24 hr later
- ▶ Complexes prepared using the LMD2 demonstrate equivalent *in vivo* gene expression to complexes prepared by hand (p>0.05 Student's t-test)

## CONCLUSIONS

- ▶ The LMD2 is a fully automated mixing device that allows reliable and reproducible production of non-viral gene therapy formulations
- ▶ Suitable for a range of different non-viral GTAs
- ▶ Easily scaleable for clinical manufacture
- ▶ Allows standardised production of non-viral gene therapy formulations for clinical applications